

BF-18M SURGI MESO THERAPY MACHINE with Pore shrinking probe

(Length: 30cm; Width:22cm; Height: 5cm & 9 cm)





Work Theory

It uses Electroporation and Electroosmosis techniques, (needle-free mesotherapy or mesoderm or cell rupture) technique, without any injectors, easily opens the aquaporin(AQP), transports skin wanted nutrition, so as to treat all kinds of skin problems. With the application of Electrophoresis energy to the beauty industry, It takes 'point to point' infiltration techniques, by orientation, layer-oriented and ration to transport high-energy nutrition serum directly into the skin tissue, let it quickly being absorbed and effective, promote skin metabolism, keep skin elastic and fine.

Function:

Freckle removal and skin whitening: Open cells channel woundless, obsorb nutrition, participate in cells metabolism, promote cells rebirth.

Fancy wrinkle removal:The electrode nutrition is directly penetrated the deep-seated skin, start cells rebirth, remove wrinkles of deep skin layer.

Resist the acne: Through the electric nutrition conveyer belt, inject the acne removal distillate into the focus, clear or restrain the acne rebirth, avoid the wound infection.

Moisten and make up water: Inject the water feeding active distillate into the deep-seated skin, supply enough water to the cells, prevent the series of skin problems caused by short of water.

ACCESSORIES

Meso Gun for nutrition/ medicine (with bottles)



Conductive stick



Pore shrinking head/ nutrition locking head



1) Mesogun:

According to the treatment of customer, to fill in relevant essence or nutrition element, they equip with the Mesogun. By operating the mesogun, to inject the essence into the deep seated skin, get the best beauty effect. The beautician can operate the Mesogun by circle gesture, move it by come and go on the treated skin surface.

2) Conductive stick:

During treatment, the probe should be touched with the customer skin fully, normally under the customer back, let the energy form a effective circuit in the human magnetic field.

3) Cooling head:

After the skin care products being injected into the skin, to avoid the nutrition releasing, can use the hammer to freeze the skin, consolidate the nutrition in the skin well.

Applicable to the crowd:

1. Face skin improving:

- (1) Oily skin, big pores, thick horniness;
- (2) Skin color dark, fatigue, pallor, dropsy and lacking oxygen.
- (3) Mature skin, slack, short of water, wrinkle, dry.

2. Eye care

- (1) Black eye socket, eye bag, wrinkle, dropsy.
- (2) Blood vessels, detoxification.